



CooperVision®  
Live Brightly.®

# Product Specifications Guide



## Proclear® family

### Proclear® multifocal



Helps provide all-day comfort for multifocal contact lens wearers who experience eye dryness.

**Wearing Schedule**  
Daily

**Carton Size Options**  
6 pack

**Material / H<sub>2</sub>O Content**  
omafilcon B / 62%

**Base Curve** 8.7mm  
**Diameter** 14.4mm

**Sphere Power**  
Plano to -8.00D  
(0.50D steps after -6.00D)  
  
+0.25D to +6.00D  
(0.25D steps)


**Add Power**  
+1.00  
+1.50  
+2.00  
+2.50

**Design**  
D lens  
N lens

**FDA Group** 2  
**UV Blocking\*** No  
**Visibility Tint** Yes

## Avaira Vitality™ family

### Avaira Vitality™



High performance. Outstanding value.

**Wearing Schedule**  
Daily

**Carton Size Options**  
3 or 6 pack

**Material / H<sub>2</sub>O Content**  
fanfilcon A / 55%


**Base Curve** 8.4mm  
**Diameter** 14.2mm

**Sphere Power**  
-0.50 to -6.00DS  
(0.25D steps)  
  
-6.50 to -12.00DS  
(0.50D steps)

**Design**  
Aspheric front surface

**FDA Group** SiHy  
**UV Blocking\*** Yes  
**Visibility Tint** Yes

### Avaira Vitality™ toric



High performance. Outstanding value.

**Wearing Schedule**  
Daily

**Carton Size Options**  
3 or 6 pack

**Material / H<sub>2</sub>O Content**  
fanfilcon A / 55%

**Base Curve** 8.5mm  
**Diameter** 14.5mm

**Sphere Power**  
+0.25 to +6.00DS  
(0.25D steps)  
  
+6.50 to +8.00DS  
(0.50D steps)  
  
-0.25 to -6.00DS  
(0.25D steps)  
  
-6.50 to -10.00DS  
(0.50D steps)

**Cylinder Power**  
-0.75  
-1.25  
-1.75  
-2.25

**Axis**  
10° to 180°  
(10° steps)

**Design**  
Optimized Toric Lens Geometry™

**FDA Group** SiHy  
**UV Blocking\*** Yes  
**Visibility Tint** Yes

## Biomedics® family

### Biomedics® Now



Designed for comfort and performance at an affordable price.

#### Wearing Schedule

Daily

#### Carton Size Options

6 pack

#### Material / H<sub>2</sub>O Content

ocufilcon D / 55%

**Base Curve** 8.8mm

**Diameter** 14.2mm

**FDA Group** 4

**UV Blocking\*** Yes

**Visibility Tint** Yes

#### Sphere Power

-0.50D to -10.00D  
(0.50D steps after  
-6.00D)

#### Design

Aspheric

\* UV-blocking contact lenses help provide protection against transmission of harmful UV radiation to eye but are not substitutes for protective UV-absorbing eyewear, such as UV absorbing goggles or sunglasses, as they do not completely cover the eye or surrounding area. Continue to use UV-absorbing eyewear as directed by your eye care professional.